

Anti-Ucp1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P25874

Item Number: CM0025358



Introduction

High-quality rabbit polyclonal antibody against human UCP1 (thermogenin), validated for Western blot, immunohistochemistry, and ELISA. Detects endogenous UCP1 in human, mouse, and rat brown adipose tissue. Ideal for metabolic studies.

[Learn More](#)

Parameter	Value
Product Name	UCP1 Rabbit pAb
Remarks/Alias	UCP; SLC25A7; UCP1
Species	Human
GeneID (Human)	7350
GeneID	7350
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human UCP1 (NP_068605.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P25874
Protein Weight	33kDa
Shipping	Ice bag