

Anti-Ucp1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If-P, Elisa - P25874

Item Number: CM0028632



Introduction

Rabbit polyclonal antibody targeting human UCP1 (Thermogenin), validated for WB, IHC-P, IF-P, and ELISA. Cross-reacts with mouse and rat. Recombinant immunogen covering aa 154-253. Ideal for mitochondrial thermogenesis and brown adipose tissue research.

[Learn More](#)

Parameter	Value
Product Name	UCP1 Rabbit pAb
Remarks/Alias	UCP; SLC25A7; UCP1
Species	Human
Gene ID	7350
Immunogen	Recombinant protein / Synthetic peptide corresponding to amino acids 154-253 of human UCP1 (NP_068605.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P25874
Protein Weight	33kDa
Shipping	Ice bag