

Anti-Ucp1+Ucp3 Rabbit Monoclonal Antibody For If-P, Ihc-P, Elisa - P12242 / P56501

Item Number: CM0003584

Introduction

Rabbit monoclonal antibody targeting mouse UCP1 and UCP3, validated for IF-P, IHC-P, and ELISA. Reacts with mouse and rat. Ideal for adipose tissue and metabolic research.

[Learn More](#)

Parameter	Value
Product Name	UCP1+UCP3 Rabbit mAb
Remarks/Alias	AI385626; Slc25a7; Ucp; UCP1+UCP3
Species	Mouse
Gene ID (Human)	7350/7352
Gene ID	22227/22229
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P12242/P56501
Protein Weight	34kDa
Shipping	Ice bag