

Anti-Ucp3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P55916

Item Number: CM0024687

Introduction

Rabbit polyclonal anti-UCP3 antibody validated for WB, IF/ICC, ELISA in human, mouse, rat. Targets the 34 kDa mitochondrial uncoupling protein 3 implicated in obesity, diabetes, and metabolic regulation.

[Learn More](#)

Parameter	Value
Product Name	UCP3 Rabbit pAb
Remarks/Alias	SLC25A9; UCP3
Species	Human
GenelD (human)	7352
GenelD	7352
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human UCP3 (NP_073714.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55916
Protein Weight	34kDa
Shipping	Ice bag