

Anti-Ue2B Monoclonal Antibody For Wb, Ihc-P, Elisa - P63146

Item Number: CM0002836



Introduction

Rabbit monoclonal antibody against human UBE2B (RAD6B), validated for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Ideal for ubiquitin pathway and DNA repair research.

[Learn More](#)

Parameter	Value
Product Name	UBE2B Rabbit mAb
Remarks/Alias	HR6B, UBC2, HHR6B, RAD6B, E2-17kDa
Species	Human
Gene ID (Human)	7320
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 53-152 of human UBE2B (P63146)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P63146
Protein Weight	17kDa
Shipping	Ice bag