

Anti-Ufc1 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Y3C8

Item Number: CM0007227



Introduction

Anti-UFC1 rabbit monoclonal antibody validated for Western blotting (WB) and ELISA.

Recognizes human, mouse, and rat UFC1, a key E2-like enzyme in the ufmylation pathway.

Suitable for research on ribosome recycling, DNA damage response, and ER-phagy.

[Learn More](#)

Parameter	Value
Product Name	UFC1 Rabbit mAb
Remarks/Alias	NEDSG; HSPC155; UFC1
Species	Human
Gene ID (Human)	51506
Gene ID	51506
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of Human UFC1 (NP_057490.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y3C8
Protein Weight	19kDa
Shipping	Ice bag