


Anti-Ufd1L Rabbit Polyclonal Antibody For Wb, Elisa - Q92890

Item Number: CM0025429



Introduction

Rabbit polyclonal antibody targeting human UFD1L (ubiquitin fusion degradation protein 1), validated for Western Blot and ELISA. Cross-reacts with mouse and rat. For ER-associated degradation and ubiquitin pathway research.

[Learn More](#)

Parameter	Value
Product Name	UFD1L Rabbit pAb
Remarks/Alias	UFD1L
Species	Human
GeneID (Human)	7353
GeneID	7353
Immunogen	Recombinant protein; Recombinant fusion protein containing amino acids 1-307 of human UFD1L (NP_005650.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92890
Protein Weight	35kDa
Shipping	Ice bag