


Anti-Ufd1L Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q92890

Item Number: CM0029689



Introduction

Rabbit polyclonal antibody targeting human UFD1L, suitable for Western blot, IF/ICC, and ELISA. Cross-reacts with human and mouse. UFD1L functions in ER-associated degradation and mitotic regulation.

[Learn More](#)

Parameter	Value
Product Name	UFD1L Rabbit pAb
Remarks/Alias	UFD1L
Species	Human
GenelD (human)	7353
GenelD	7353
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q92890
Protein Weight	35kDa
Shipping	Ice bag