


Anti-Ufl1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O94874

Item Number: CM0021580



Introduction

Rabbit polyclonal antibody targeting human UFL1 (E3 UFM1-protein ligase 1). Validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for ufmylation, ER stress and DNA damage research.

[Learn More](#)

Parameter	Value
Product Name	UFL1 Rabbit pAb
Remarks/Alias	NLBP; RCAD; Maxer; KIAA0776; UFL1
Species	Human
GeneID (Human)	23376
GeneID	23376
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 601-794 of human UFL1 (NP_056138.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O94874
Protein Weight	90kDa
Shipping	Ice bag