


Anti-Ufm1 Monoclonal Antibody For Wb, Ihc-P, Elisa - P61960

Item Number: CM0002216



Introduction

Rabbit monoclonal antibody targeting human UFM1 (ubiquitin-fold modifier 1). Validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat. Ideal for studying ufmylation in DNA damage, ER stress, and ribosome recycling.

[Learn More](#)

Parameter	Value
Product Name	UFM1 Rabbit mAb
Remarks/Alias	HLD14; BM-002; C13orf20; UFM1
Species	Human
Gene ID (Human)	51569
Gene ID	51569
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-85 of human UFM1 (P61960).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61960
Protein Weight	9kDa
Shipping	Ice bag