


Anti-Ufm1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P61960

Item Number: CM0021867



Introduction

Rabbit polyclonal antibody targeting human UFM1, validated in WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Used for ufmylation, ribosome biology, and ER stress research.

[Learn More](#)

Parameter	Value
Product Name	UFM1 Rabbit pAb
Remarks/Alias	HLD14; BM-002; C13orf20; UFM1
Species	Human
Gene ID (Human)	51569
Gene ID	51569
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-63 of human UFM1 (NP_057701.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61960
Protein Weight	9kDa
Shipping	Ice bag