


## Anti-Ugdh Polyclonal Antibody For Ihc-P And Elisa - O60701

Item Number: CM0018574



### Introduction

Rabbit polyclonal antibody targeting UGDH (UDP-glucose 6-dehydrogenase), validated for immunohistochemistry (IHC-P) and ELISA. Reacts with human, mouse, and rat samples. Useful for glycosaminoglycan biosynthesis and neurodevelopmental disease research.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | UGDH Rabbit pAb                                                                                                   |
| Remarks/Alias    | GDH; UGD; DEE84; EIEE84; UDPGDH; UDP-GlcDH; UGDH                                                                  |
| Species          | Human                                                                                                             |
| GeneID (human)   | 7358                                                                                                              |
| GeneID           | 7358                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 165-494 of human UGDH (NP_003350.1) |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | IHC-P, ELISA                                                                                                      |
| Cross Reactivity | Human, Mouse, Rat                                                                                                 |
| SWISS            | O60701                                                                                                            |
| Protein Weight   | 55kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |