


Anti-Ugdh Rabbit Polyclonal Antibody For Wb And Elisa - O60701

Item Number: CM0014653



Introduction

Rabbit polyclonal antibody against UGDH (UDP-glucose 6-dehydrogenase), validated for Western blot and ELISA. Detects human and mouse UGDH, a key enzyme in glycosaminoglycan synthesis, with a predicted molecular weight of 55 kDa.

[Learn More](#)

Parameter	Value
Product Name	UGDH Rabbit pAb
Remarks/Alias	GDH; UGD; DEE84; EIEE84; UDPGDH; UDP-GlcDH; UGDH
Species	Human
GeneID (Human)	7358
GeneID	7358
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O60701
Protein Weight	55kDa
Shipping	Ice bag