

Anti-Uggt1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nyu2

Item Number: CM0028176

Introduction

High-quality Anti-UGGT1 Rabbit Polyclonal Antibody validated for WB and ELISA. Detects human, mouse, rat UGGT1 (177 kDa), an ER glycoprotein folding sensor. Ideal for ER stress and protein quality control research.

[Learn More](#)

Parameter	Value
Product Name	UGGT1 Rabbit pAb
Remarks/Alias	UGT1; HUGT1; UGCGL1; UGGT1
Species	Human
Gene ID (Human)	56886
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1456-1555 of human UGGT1 (NP_064505.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NYU2
Protein Weight	177kDa
Shipping	Ice bag