

Anti-Ugt1A1 Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P22309

Item Number: CM0029039

Introduction

This rabbit polyclonal antibody against human UGT1A1 (SWISS: P22309) is suitable for WB, IF/ICC, IHC-P, and ELISA. Validated in human, mouse, rat. Ideal for UDP-glucuronosyltransferase research.

[Learn More](#)

Parameter	Value
Product Name	UGT1A1 Rabbit pAb
Remarks/Alias	GNT1; UGT1; UDPGT; UGT1A; HUG-BR1; BILIQTL1; UDPGT 1-1; UGT1A1
Species	Human
Gene ID (Human)	54658
Gene ID	54658
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human UGT1A1 (NP_000454.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22309
Protein Weight	60kDa
Shipping	Ice bag