

Anti-Ugt1A1 Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P22309

Item Number: CM0017770

Introduction

Rabbit polyclonal UGT1A1 antibody validated for WB, IF/ICC, IHC-P, and ELISA. Detects UDP-glucuronosyltransferase 1A1 in human, mouse, and rat samples. Ideal for drug metabolism, bilirubin conjugation, and toxicity studies.

[Learn More](#)

Parameter	Value
Product Name	UGT1A1 Rabbit pAb
Remarks/Alias	GNT1; UGT1; UDPGT; UGT1A; HUG-BR1; BILIQL1; UDPGT 1-1; UGT1A1
Species	Human
GeneID (Human)	54658
GeneID	54658
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human UGT1A1 (NP_000454.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22309
Protein Weight	60 kDa
Shipping	Ice bag