

Anti-Ugt1A1 Rabbit Monoclonal Antibody For Wb, Elisa - P22309

Item Number: CM0006840

Introduction

Rabbit monoclonal antibody against human UGT1A1, suitable for Western blot and ELISA, cross-reacts with rat. Relevant for phase II drug metabolism, bilirubin conjugation, and hormone glucuronidation studies.

[Learn More](#)

Parameter	Value
Product Name	UGT1A1 Rabbit mAb
Remarks/Alias	GNT1; UGT1; UDPGT; UGT1A; HUG-BR1; BILIQTL1; UDPGT 1-1; UGT1A1
Species	Human
Gene ID (Human)	54658
Gene ID	54658
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P22309
Protein Weight	60 kDa
Shipping	Ice bag