

Anti-Ugt1A6 Polyclonal Antibody For Wb, Elisa - P19224

Item Number: CM0023308

Introduction

Rabbit polyclonal antibody against human UGT1A6 (UniProt P19224) for Western blot and ELISA applications. Detects endogenous UGT1A6 in human and mouse samples. Suitable for drug metabolism and detoxification research.

[Learn More](#)

Parameter	Value
Product Name	UGT1A6 Rabbit pAb
Remarks/Alias	GNT1, UGT1, HLUGP, UDPGT, UGT1A, UGT1C, UGT1E, UGT1F, HLUGP1, UGT-1A, UGT-1C, UGT-1E, UGT-1F, UGT1.1, UGT1.3, UGT1.5, UGT1.6, UGT1A1, UGT1A3, UGT1A5, UGT1-01, UGT1-03, UGT1-05, UGT1-06, UGT1A6S, hUG-BR1, UDPGT 1-6, UGT1A6
Species	Human
Gene ID (Human)	54578
Immunogen	Recombinant Protein corresponding to a sequence within amino acids 65-270 of human UGT1A6 (NP_001063.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P19224
Protein Weight	61kDa
Shipping	Ice bag