

Anti-Ugt1A6 Polyclonal Antibody For Wb, If-P, Elisa - P19224

Item Number: CM0020102

Introduction

Rabbit polyclonal antibody against human UGT1A6, validated for Western blot, immunofluorescence, and ELISA. Suitable for human, mouse, and rat samples. Ideal for drug metabolism and detoxification research.

[Learn More](#)

Parameter	Value
Product Name	UGT1A6 Rabbit pAb
Remarks/Alias	GNT1; UGT1; HLUGP; UDPGT; UGT1A; UGT1C; UGT1E; UGT1F; HLUGP1; UGT-1A; UGT-1C; UGT-1E; UGT-1F; UGT1.1; UGT1.3; UGT1.5; UGT1.6; UGT1A1; UGT1A3; UGT1A5; UGT1-01; UGT1-03; UGT1-05; UGT1-06; UGT1A6S; HUG-BR1; UDPGT 1-6; UGT1A6
Species	Human
GeneID (Human)	54578
GeneID	54578
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 65-270 of human UGT1A6 (NP_001063.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19224
Protein Weight	61kDa
Shipping	Ice bag