

Anti-Ugt2B10 Rabbit Polyclonal Antibody For Wb And Elisa - P36537

Item Number: CM0028693



Introduction

Rabbit anti-UGT2B10 polyclonal antibody validated for WB and ELISA. Recognizes human, mouse, and rat UGT2B10 (61 kDa), a key phase II UDP-glucuronosyltransferase. Ideal for drug metabolism and detoxification research.

[Learn More](#)

Parameter	Value
Product Name	UGT2B10 Rabbit pAb
Remarks/Alias	UDPGT2B10; UGT2B10
Species	Human
GeneID (Human)	7365
GeneID	7365
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P36537
Protein Weight	61kDa
Shipping	Ice bag