


Anti-Ugt2B17 Rabbit Polyclonal Antibody For If/Icc And Elisa - O75795

Item Number: CM0019624



Introduction

Rabbit polyclonal antibody targeting human UGT2B17, a key enzyme in steroid hormone glucuronidation. Validated for IF/ICC and ELISA. Suitable for drug metabolism and endocrine research studies.

[Learn More](#)

Parameter	Value
Product Name	UGT2B17 Rabbit pAb
Remarks/Alias	BMND12; UDPGT2B17; UGT2B17
Species	Human
GeneID (Human)	7367
GeneID	7367
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-200 of human UGT2B17 (NP_001068.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	O75795
Protein Weight	61kDa
Shipping	Ice bag