

Anti-Ugt3A2 Polyclonal Antibody For Wb, Elisa - Q3Sy77

Item Number: CM0030780

Introduction

Rabbit polyclonal antibody against human UGT3A2, a UDP-glycosyltransferase involved in xenobiotic metabolism. Suitable for WB and ELISA, with cross-reactivity to mouse and rat. UniProt: Q3SY77.

[Learn More](#)

Parameter	Description
Product Name	UGT3A2 Rabbit pAb
Remarks/Alias	UGT3A2
Species	Human
Gene ID	167127
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-300 of human UGT3A2 (NP_777574.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q3SY77
Protein Weight	60kDa
Shipping	Ice bag