

Anti-Uhrf1 Polyclonal Antibody For Wb And Elisa - Q96T88

Item Number: CM0030537



Introduction

Rabbit polyclonal antibody targeting human UHRF1 (ICBP90/Np95), validated for western blot and ELISA. Ideal for epigenetics research, including DNA methylation, histone modifications, and cancer biology.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] UHRF1 Rabbit pAb
Remarks/Alias	Np95; hNP95; ICBP90; RNF106; TDRD22; hUHRF1; huNp95
Species	Human
GeneID (Human)	29128
GeneID	29128
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 672-793 of human UHRF1 (NP_037414.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96T88
Protein Weight	90kDa
Shipping	Ice bag