

Anti-Uhrf1 Polyclonal Antibody For Wb, Elisa - Q96T88

Item Number: CM0022604



Introduction

High-quality UHRF1 rabbit polyclonal antibody for Western blot and ELISA applications, recognizing human and mouse targets. Essential for DNA methylation and chromatin modification research.

[Learn More](#)

Parameter	Value
Product Name	UHRF1 Rabbit pAb
Remarks/Alias	Np95; hNP95; ICBP90; RNF106; TDRD22; hUHRF1; huNp95; UHRF1
Species	Human
Gene ID	29128 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human UHRF1 (NP_001041666.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96T88
Protein Weight	90kDa
Shipping	Ice bag