

Anti-Uhrf1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q96T88

Item Number: CM0011823



Introduction

Rabbit monoclonal antibody against human UHRF1 (SWISS Q96T88). Reacts with human, mouse, rat. Validated for WB, IHC-P, ELISA. Key epigenetic regulator in DNA methylation inheritance and chromatin modification.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] UHRF1 Rabbit mAb
Remarks/Alias	Np95; hNP95; ICBP90; RNF106; TDRD22; hUHRF1; huNp95
Species	Human
Gene ID (Human)	29128
Gene ID	29128
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human UHRF1 (NP_001276981.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96T88
Protein Weight	90kDa
Shipping	Ice bag