

Anti-Uhrf2 Rabbit Polyclonal Antibody For Wb And Elisa - Q7Tmi3

Item Number: CM0023127

Introduction

High-quality rabbit polyclonal antibody against mouse UHRF2, validated for Western blot and ELISA. Detects human, mouse, rat UHRF2, a key E3 ubiquitin ligase in DNA methylation and cell cycle regulation. Ideal for epigenetics and DNA repair research.

[Learn More](#)

Parameter	Value
Product Name	Uhrf2 Rabbit pAb
Remarks/Alias	Nirf; 2310065A22Rik; D130071B19Rik; Uhrf2
Species	Mouse
GeneID (Human)	115426
GeneID	109113
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 76-359 of mouse Uhrf2 (NP_659122.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7TMI3
Protein Weight	90kDa
Shipping	Ice bag