

Anti-Uimc1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q96RI1

Item Number: CM0028450



Introduction

Rabbit polyclonal antibody targeting human UIMC1 (RAP80), validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Immunogen: N-terminal 1-150 aa. Suitable for DNA damage and repair research.

[Learn More](#)

Parameter	Value
Product Name	UIMC1 Rabbit pAb
Remarks/Alias	RAP80; X2HRIP110; UIMC1
Species	Human
GeneID (Human)	51720
GeneID	51720
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human UIMC1 (NP_057374.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96RL1
Protein Weight	80kDa
Shipping	Ice bag