

Anti-Ulbp1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Bzm6

Item Number: CM0016731

Introduction

High-quality rabbit polyclonal antibody against human ULBP1 (NKG2D ligand 1). Validated for WB, IHC-P, ELISA. Cross-reacts with mouse. Ideal for immune surveillance and NK cell research.

[Learn More](#)

Parameter	Specification
Product Name	ULBP1 Rabbit pAb
Remarks/Alias	N2DL-1; RAET1I; NKG2DL1; ULBP1
Species	Human
Gene ID	80329 (Human)
Immunogen	Recombinant fusion protein containing amino acids 26-216 of human ULBP1 (NP_079494.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BZM6
Protein Weight	28kDa
Shipping	Ice bag