

Anti-ULbp2 Rabbit Monoclonal Antibody For Wb And Elisa - Q9Bzm5

Item Number: CM0007019



Introduction

Rabbit monoclonal antibody against human ULBP2 (NKG2D ligand), validated for Western blot and ELISA. Reacts specifically with human ULBP2, a 27 kDa GPI-anchored protein that activates NK cell cytotoxicity. Ideal for cancer and immunology research.

[Learn More](#)

Parameter	Value
Product Name	ULBP2 Rabbit mAb
Remarks/Alias	ULBP2; ALCAN-alpha; N2DL2; NKG2DL2; RAET1H; NKG2D ligand 2
Species	Human
GeneID (Human)	80328
GeneID	80328
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-217 of human ULBP2 (NP_079493.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9BZM5
Protein Weight	27kDa
Shipping	Ice bag