

Anti-ULK1 (Phospho S555) Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - O75385

Item Number: CM0001264



Introduction

Phospho-ULK1-S555 Rabbit mAb specifically detects ULK1 phosphorylated at Ser555. Validated for WB, IF/ICC, and ELISA. Reacts with human and mouse. Ideal for autophagy and signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-ULK1-S555 Rabbit mAb
Remarks/Alias	ATG1; ATG1A; UNC51; hATG1; Unc51.1
Species	Human
GeneID (Human)	8408
GeneID	8408
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S555 of human ULK1 (NP_003556.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O75385
Protein Weight	113kDa
Shipping	Ice bag