

Anti-ULK1 Rabbit Polyclonal Antibody For Wb And Elisa - O75385

Item Number: CM0019094



Introduction

Rabbit polyclonal antibody against human ULK1 (ATG1), validated in WB and ELISA, cross-reacts with mouse. For autophagy and neuronal differentiation research.

[Learn More](#)

Parameter	Value
Product Name	ULK1 Rabbit pAb
Remarks/Alias	ATG1; ATG1A; UNC51; hATG1; Unc51.1; ULK1
Species	Human
GeneID (Human)	8408
GeneID	8408
Immunogen	Recombinant protein of human ULK1.
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	O75385
Protein Weight	113kDa
Shipping	Ice bag