

Anti-ULK3 Monoclonal Antibody For Wb, If/Icc, Elisa - Q6Phr2

Item Number: CM0012073



Introduction

Recombinant rabbit monoclonal antibody against human ULK3, validated for WB, IF/ICC, and ELISA. Targets the serine/threonine kinase involved in hedgehog signaling and autophagy. Cross-reacts with human and mouse.

[Learn More](#)

Parameter	Value
Product Name	ULK3 Rabbit mAb
Remarks/Alias	ULK3
Species	Human
GeneID (Human)	25989
GeneID	25989
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 373-472 of human ULK3 (Q6PHR2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q6PHR2
Protein Weight	53kDa
Shipping	Ice bag