

Anti-Umod Rabbit Polyclonal Antibody For Wb, Elisa - P07911

Item Number: CM0023929



Introduction

UMOD Rabbit Polyclonal Antibody validated for WB and ELISA with reactivity in Human and Mouse. Ideal for studying uromodulin function in kidney apical membrane biology and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	UMOD Rabbit pAb
Remarks/Alias	THP; FJHN; HNFJ; THGP; HNFJ1; MCKD2; ADTKD1; ADMCKD2; UMOD
Species	Human
GeneID (Human)	7369
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P07911
Protein Weight	70kDa
Shipping	Ice bag