

## Anti-Umps Polyclonal Antibody For Wb, Elisa - P11172

Item Number: CM0025259



### Introduction

Rabbit polyclonal antibody targeting human Umps (OPRT/ODC). Validated for WB and ELISA. Reacts with human samples. Detects a band at 52 kDa. Ideal for pyrimidine biosynthesis research and IVD assay development.

[Learn More](#)

| Parameter        | Value                                                                                                                                  |
|------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | UMPS Rabbit pAb                                                                                                                        |
| Remarks/Alias    | OPRT; UMPS                                                                                                                             |
| Species          | Human                                                                                                                                  |
| GeneID (human)   | 7372                                                                                                                                   |
| GeneID           | 7372                                                                                                                                   |
| Immunogen        | Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 181-480 of human UMPS (NP_000364.1) |
| Source           | Rabbit                                                                                                                                 |
| Category         | Polyclonal Antibodies                                                                                                                  |
| Application      | WB, ELISA                                                                                                                              |
| Cross Reactivity | Human                                                                                                                                  |
| SWISS            | P11172                                                                                                                                 |
| Protein Weight   | 52kDa                                                                                                                                  |
| Shipping         | Ice bag                                                                                                                                |