

## Anti-Umps Polyclonal Antibody For Wb, If/Icc, Elisa - P11172

Item Number: CM0017564



### Introduction

Rabbit polyclonal antibody targeting human UMPS, validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for pyrimidine biosynthesis research.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | UMPS Rabbit pAb                                                                                                    |
| Remarks/Alias    | OPRT; UMPS                                                                                                         |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 7372                                                                                                               |
| GeneID           | 7372                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 100-260 of human UMPS (NP_000364.1). |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, IF/ICC, ELISA                                                                                                  |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | P11172                                                                                                             |
| Protein Weight   | 52kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |