

Anti-Unc13B Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O14795

Item Number: CM0021865

Introduction

High-quality rabbit anti-UNC13B polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human and mouse UNC13B, a synaptic vesicle priming protein involved in exocytosis. Suitable for neuroscience and cell biology research.

[Learn More](#)

Parameter	Value
Product Name	UNC13B Rabbit pAb
Remarks/Alias	UNC13; MUNC13; Unc13h2; munc13-2; UNC13B
Species	Human
Gene ID (Human)	10497
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-350 of human UNC13B (NP_006368.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O14795
Protein Weight	181kDa
Shipping	Ice bag