

Anti-Unc13D Polyclonal Antibody For Wb, Elisa - Q70J99

Item Number: CM0014693



Introduction

Rabbit polyclonal antibody against human UNC13D (Munc13-4) for WB and ELISA, cross-reactive with mouse. Validated for research in cytotoxic granule exocytosis and immune cell function.

[Learn More](#)

Parameter	Value
Product Name	UNC13D Rabbit pAb
Remarks/Alias	FHL3, HLH3, HPLH3, Munc13-4, UNC13D
Species	Human
GeneID (Human)	201294
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 891-1090 of human UNC13D (NP_954712.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q70J99
Protein Weight	123kDa
Shipping	Ice bag