

Anti-Unc50 Rabbit Polyclonal Antibody For Wb, Elisa - Q53Hi1

Item Number: CM0020441

Introduction

Rabbit polyclonal antibody against human UNC50, validated for WB and ELISA, reacts with human, mouse, and rat. UNC50 is involved in neuronal nicotinic receptor expression and RNA binding.

[Learn More](#)

Parameter	Value
Product Name	UNC50 Rabbit pAb
Remarks/Alias	URP; GMH1; UNCL; HSD23; PDLs22; UNC50
Species	Human
GeneID (Human)	25972
GeneID	25972
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q53HI1
Protein Weight	30kDa
Shipping	Ice bag