

Anti-Unc5B Polyclonal Antibody For Wb, Elisa - Q8Izj1

Item Number: CM0018316

Introduction

Validated UNC5B rabbit polyclonal antibody for WB and ELISA, reactive with human, mouse, rat. Ideal for netrin receptor research in axon guidance, angiogenesis, and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	UNC5B Rabbit pAb
Remarks/Alias	UNC5H2; p53RDL1; UNC5B
Species	Human
GeneID (Human)	219699
GeneID	219699
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 290-377 of human UNC5B (NP_734465.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IZJ1
Protein Weight	104kDa
Shipping	Ice bag