

Anti-Unc93B1 Polyclonal Antibody For Wb, Elisa - Q9H1C4

Item Number: CM0019945



Introduction

High-quality Anti-UNC93B1 Rabbit Polyclonal Antibody validated for WB and ELISA. Suitable for innate immunity and TLR signaling research. Reacts with Human, Mouse, Rat. Ideal for diagnostic and research labs.

[Learn More](#)

Parameter	Value
Product Name	UNC93B1 Rabbit pAb
Remarks / Alias	IIAE1; UNC93; UNC93B; Unc-93B1; UNC93B1
Species	Human
Gene ID (Human)	81622
Gene ID	81622
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 450-550 of human UNC93B1 (NP_112192.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H1C4
Protein Weight	67kDa
Shipping	Ice bag