

Anti-Upb1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Ubr1

Item Number: CM0021710



Introduction

Rabbit polyclonal antibody targeting human UPB1 (beta-ureidopropionase), validated for Western blot and ELISA. Cross-reacts with mouse and rat. Suitable for pyrimidine degradation and metabolic disorder research.

[Learn More](#)

Parameter	Value
Product Name	UPB1 Rabbit pAb
Remarks/Alias	BUP-1; UPB1
Species	Human
GeneID (Human)	51733
GeneID	51733
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human UPB1 (NP_057411.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9UBR1
Protein Weight	43kDa
Shipping	Ice bag