

Anti-Upf1/Rent1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q92900

Item Number: CM0015058



Introduction

Rabbit polyclonal antibody against UPF1/RENT1, validated for WB, IHC-P, ELISA. Reacts with human and mouse. Ideal for nonsense-mediated decay research. Immunogen: synthetic peptide. ~124 kDa.

[Learn More](#)

Parameter	Value
Product Name	UPF1/RENT1 Rabbit pAb
Remarks/Alias	UTF; HUPF1; NORF1; RENT1; smg-2; pNORF1; UPF1/RENT1
Species	Human
GeneID (Human)	5976
GeneID	5976
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q92900
Protein Weight	124kDa
Shipping	Ice bag