

Anti-Upf1/Rent1 Rabbit Monoclonal Antibody For Wb, If/Icc, Ip, Elisa - Q92900

Item Number: CM0011449



Introduction

Rabbit monoclonal antibody targeting human UPF1/RENT1, validated for WB, IF/ICC, IP, and ELISA. Cross-reacts with mouse and rat. Suitable for nonsense-mediated mRNA decay and RNA helicase studies. Molecular weight: 124 kDa.

[Learn More](#)

Parameter	Value
Product Name	UPF1/RENT1 Rabbit mAb
Remarks/Alias	UTF; HUPF1; NORF1; RENT1; smg-2; pNORF1; UPF1/RENT1
Species	Human
Gene ID (Human)	5976
Gene ID	5976
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human UPF1/RENT1 (Q92900)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92900
Protein Weight	124kDa
Shipping	Ice bag