


Anti-Upf2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q9Hau5

Item Number: CM0030644



Introduction

Rabbit polyclonal antibody against UPF2, a key regulator of nonsense-mediated mRNA decay. Validated for WB, IHC-P, IP, ELISA. Reacts with human and mouse. Ideal for NMD research.

[Learn More](#)

Parameter	Value
Product Name	UPF2 Rabbit pAb
Remarks/Alias	HUPF2; RENT2; smg-3; UPF2
Species	Human
Gene ID (Human)	26019
Gene ID	26019
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1153-1272 of human UPF2 (NP_542166.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9HAU5
Protein Weight	148kDa
Shipping	Ice bag