


Anti-Upf2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q9Hau5

Item Number: CM0014498



Introduction

High-quality rabbit polyclonal antibody targeting human UPF2 (Q9HAU5), validated in WB, IHC-P, IP, and ELISA. Ideal for nonsense-mediated mRNA decay (NMD) research. Cross-reacts with human and rat.

[Learn More](#)

Parameter	Value
Product Name	UPF2 Rabbit pAb
Remarks/Alias	HUPF2; RENT2; smg-3; UPF2
Species	Human
GeneID (Human)	26019
GeneID	26019
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1153-1272 of human UPF2 (NP_542166.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9HAU5
Protein Weight	148kDa
Shipping	Ice bag