

Anti-Upk1B Polyclonal Antibody For Western Blot, Elisa - O75841

Item Number: CM0026963

Introduction

Anti-UPK1B rabbit polyclonal antibody validated for Western blot and ELISA. Recognizes human and mouse uroplakin-1b (UPIB/TSPAN20), suitable for bladder epithelial research and urothelial plaque studies.

[Learn More](#)

Parameter	Value
Product Name	UPK1B Rabbit pAb
Remarks/Alias	UPIB; UPK1; TSPAN20; UPK1B
Species	Human
GeneID (Human)	7348
GeneID	7348
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 110-230 of human UPK1B (NP_008883.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O75841
Protein Weight	30kDa
Shipping	Ice bag