

Anti-Upk2 Rabbit Polyclonal Antibody For Wb, Elisa - 000526

Item Number: CM0015412

Introduction

High-quality rabbit polyclonal antibody against human Uroplakin-2 (UPK2) validated for WB and ELISA. Recombinant protein immunogen (aa 85-155) ensures specific detection. Ideal for urothelial asymmetric unit membrane research.

[Learn More](#)

Parameter	Value
Product Name	UPK2 Rabbit pAb
Remarks/Alias	UP2; UPII; UPK2
Species	Human
GeneID (Human)	7379
GeneID	7379
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 85-155 of human UPK2 (NP_006751.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O00526
Protein Weight	19kDa
Shipping	Ice bag