

Anti-Upk3A Rabbit Polyclonal Antibody For Wb, If-P, Elisa - O75631

Item Number: CM0013397

Introduction

High-quality rabbit polyclonal antibody against human UPK3A (uroplakin-3a). Validated for WB, IF-P, ELISA. Reacts with human, mouse, rat. Ideal for urothelial cell differentiation studies.

[Learn More](#)

Parameter	Value
Product Name	UPK3A Rabbit pAb
Remarks/Alias	UP3A; UPK3; UPIII; UPIIIA; UPK3A
Species	Human
Gene ID (Human)	7380
Gene ID	7380
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-207 of human UPK3A (NP_008884.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75631
Protein Weight	31kDa
Shipping	Ice bag