

Anti-Upk3B Rabbit Polyclonal Antibody For Wb, Elisa - Q9Bt76

Item Number: CM0020896

Introduction

UPK3B Rabbit pAb is a polyclonal antibody raised in rabbit, validated for WB and ELISA in human samples, targeting uroplakin-3b (Q9BT76). Ideal for AUM and urothelial differentiation research.

[Learn More](#)

Parameter	Value
Product Name	UPK3B Rabbit pAb
Remarks/Alias	P35; UP3B; UPIIB; UPK3B
Species	Human
GeneID (Human)	105375355
GeneID	105375355
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 60-200 of human UPK3B (NP_085047.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9BT76
Protein Weight	34kDa
Shipping	Ice bag