

Anti-Upp1 Rabbit Polyclonal Antibody For Wb And Elisa - Q16831

Item Number: CM0015534



Introduction

Rabbit polyclonal antibody targeting human uridine phosphorylase 1 (UPP1). Suitable for Western blot and ELISA; detects a 34 kDa protein involved in pyrimidine metabolism and nucleotide synthesis. Ideal for cancer and metabolism research.

[Learn More](#)

Parameter	Value
Product Name	UPP1 Rabbit pAb
Remarks/Alias	UP; UPP; UPASE; UDRPASE; UPP1
Species	Human
Gene ID (Human)	7378
Gene ID	7378
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 89-173 of human UPP1 (NP_001274357.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q16831
Protein Weight	34kDa
Shipping	Ice bag